

# COMPANY NEWSLETTER

[WWW.MEHROTRABIOTECH.COM](http://WWW.MEHROTRABIOTECH.COM)



September, 2023

Visit us at

Hall no. 2 | Booth no. P-25

14-16 September, 2023  
HITEX, Hyderabad



analytica Anacon India



INDIA LAB EXPO  
2023

International trade fair for Laboratory Technology,  
Analysis, Biotechnology and Diagnostics

#BiggerThanEver

[analyticalindia.com](http://analyticalindia.com) | [www.indialabexpo.com](http://www.indialabexpo.com)

## Top News

New Product Launch: Minus 80 Degree Celsius Ultra Low Temperature Deep Freezer

Mehrotra Biotech Private Limited has been nominated for 2023 Go Global Awards

We are participating in INDIA LAB EXPO-  
analytica Anacon India-2023

Special Features about -20 Deg C  
Deep Freezer

"UDAAN"--Annual Conference Meet  
2023-24

## Special features of minus 80 deg C Ultra Low Temperature Deep Freezer

*We are currently a well regarded manufacturer of -20 Degree Celsius Deep freezers. Moving forward, we are pleased to announce the introduction of a new line of -80 DEGREE CELSIUS ULTRA LOW TEMPERATURE DEEP FREEZER.*

The new designed assisting handle and vacuum release port let the freezer be operated singlehanded to open and close the outer door.

The 6 sides of the cabinet are made from high-performance VIP: Vacuum Insulation Panel material, improving thermal insulation performance to a large extent.

Data can be recorded every 5 minutes for up to 10 years.



The perfect audible & visual alarm system : high and low temperature ,sensor failure ,power failure ,low battery ,door ajar ,main board communication error, high ambient temperature etc.

Optional Accessories:  
WIFI/Bluetooth, Printer/Mechanical, Lock/Swipe/Card/Fingerprint

01/02

Advance Technologies  
Effective Solutions.

# COMPANY NEWSLETTER

WWW.MEHROTRABIOTECH.COM



September, 2023

MEHROTRA  
Biotech

MEHROTRA  
BIOTECH PVT LTD  
has been  
nominated for the  
2023 Go Global  
Awards

Entries Close 30th June 2023

[www.goglobalawards.org](http://www.goglobalawards.org)



## -20 Deg C Deep Freezer

Available with Glass Door &  
Solid Door with different capacities



- Precise temperature control
- High storage capacity
- Efficient insulation
- Advanced cooling system
- Alarm and monitoring features
- User-friendly interface
- Long-term storage capabilities
- Easy maintenance
- PC monitoring compatibility
- Data logging and retrieval



Mehrotra Biotech Private Limited Annual Conference was held in Goa from 14th-18th June with the theme of "UDAAN" in which a constructive conversation was held on achievements, future goals and endeavors.



02/02

Advance Technologies  
Effective Solutions.